

# Cavitation Magic Pink - User Manual & Technical Catalog

**Welcome to the Comprehensive Guide for the Cavitation Magic Pink Machine.**  
This premium body contouring and skin rejuvenation multi-functional device combines five advanced clinical technologies to deliver non-invasive fat reduction, skin tightening, and muscle stimulation.

## 1. Key Integrated Technologies

| Technology           | Primary Function & Description                                                      |
|----------------------|-------------------------------------------------------------------------------------|
| Cavitation           | Ultrasonic 40KHz waves targeting and breaking down stubborn fat cells.              |
| RF (Radio Frequency) | Thermal energy deep into dermal layers to stimulate collagen and tighten skin.      |
| Vacuum               | Suction therapy to enhance lymphatic drainage and reduce cellulite appearance.      |
| Laser (Lipo Laser)   | Low-level cold laser photobiomodulation to liquefy fat cell contents safely.        |
| Current (EMS)        | Electrical Muscle Stimulation to tone, firm, and lift underlying muscle structures. |

## 2. General Clinical Guidelines & Protocols

- **Hydration:** Ensure patients drink at least 1.5 to 2 liters of water before and after sessions to facilitate fat clearance through the lymphatic system.
- **Session Duration:** Do not exceed 30 minutes per area or a total session time of 60 minutes.
- **Conductive Media:** Always apply high-quality ultrasound gel for cavitation, and specialized radiofrequency cream or oil for RF handles to protect the epidermis.

This is for informational purposes only. For medical advice or diagnosis, consult a professional. AI responses may include mistakes.